


pLenti-GFP-Puro+IGF1R(h)质粒

Cat No.: VGP-3480-34F8

Gene/insert name	Homo sapiens insulin like growth factor 1 receptor (IGF1R)
Alt name	CD221; IGFIR; IGFR; JTK13
Insert size	4104 bp
Gene ID	3480
Accessions	NM_000875
Tag	W/O
Backbone size	8190 bp
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	GFP

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
Note	For research use only, not for use in diagnostic procedure.

