


**pLenti-GFP-Puro+shRNA-USP22(h) 五靶点干扰质粒套餐**

Cat No.: VGP-23326-SHK2B

<b>Target gene name</b>	Homo sapiens ubiquitin specific peptidase 22
<b>Official Symbol</b>	USP22
<b>Alt name</b>	USP3L
<b>Gene ID</b>	<a href="#">23326</a>
<b>Plasmid size</b>	7970 bp
<b>Cloning site 5'</b>	Age I
<b>Cloning site 3'</b>	EcoR I
<b>Promoter</b>	U6
<b>Loop</b>	CTCGAG
<b>Growth strain(s):</b>	Stbl 3
<b>High or low copy</b>	High
<b>Bacterial resistance(s)</b>	Ampicillin
<b>Growth temperature</b>	37 °C
<b>Vector type</b>	Lentiviral, Mammalian Expression
<b>Selectable markers</b>	Puromycin
<b>Fluorescent protein</b>	GFP

**Map**



<b>Buffer</b>	TE buffer
<b>Storage</b>	Store at -20°C up to 24 months
<b>The plasmid is ready for</b>	Transfection
<b>Note</b>	For research use only, not for use in diagnostic procedure.



---

Vigen Biotechnology (Zhenjiang) co., Ltd  
Technical Support: 183 628 99236  
[www.vigenbio.com](http://www.vigenbio.com)

