


pLenti-mCherry-BSD+FAP (h) plasmid

Cat No.: VGP-2191-33B8

Gene/insert name	Homo sapiens fibroblast activation protein alpha
Official Symbol	FAP
Alt name	DPPIV; FAPA; FAPalpha; SIMP
Insert size	2280 bp
Gene ID	2191
Accessions	NM_004460
Tag	W/O
Backbone size	7949 bp
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Blasticidin
Fluorescent protein	mCherry

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

