


pLenti-Puro-shRNA·NUP85(h5)靶点质粒套餐

Cat No.: VGP-79902-SHK2

Target gene	Homo sapiens nucleoporin 85
Official symbol	NUP85
Alt name	NUP85 , FROUNT; NPHS17; Nup75
Gene ID	<u>79902</u>
Tag	W/O
Backbone size	8155 bp
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	GFP
Promoter	U6
Loop	CTCGAG

Backbone Map



The plasmid is ready for	Transfection
Storage	Store at -80°C up to 12 months
Note	For research use only, not for use in diagnostic procedure.

