


pLenti-GFP-Puro+EGFRvIII^h plasmid

Cat No.: VGP-1956vIII-33C8

Gene/insert name	Homo sapiens epidermal growth factor receptor (EGFR)
Official Symbol	EGFRvIII ^h
Alt name	ERBB; ERBB1; ERBB
Insert size	1203 bp
Gene ID	1956
Accessions	NM_001346941.2
Tag	W/O
Backbone size	8194 bp
Cloning site 5'	Xho I
Cloning site 3'	EcoR I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	GFP

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

