


pLenti-Puro+mCCR5·APEX2·FLAG plasmid

Cat No.: VGP-12774-33B3

Gene/insert name	Mus musculus C-C motif chemokine receptor 5
Official Symbol	Ccr5
Alt name	AM4-7; CD195; Cmkbr5
Insert size	1062 bp
Gene ID	12774
Accessions	NM_009917.5
Tag	FLAG
Backbone size	8174 bp
Cloning site 5'	Xba I
Cloning site 3'	Xho I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	W/O

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

