


pLenti-BSD+Cas9 plasmid

Cat No.: VGP-901176-34E8

Gene/insert name	Cas9
Official Symbol	Cas9
Alt name	N/A
Insert size	4104 bp
Gene ID	N/A
Accessions	N/A
Tag	P2A
Backbone size	8651 bp
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Blasticidin
Fluorescent protein	W/O

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

