


pLenti-CMV-MCS-EF1a-mCherry-Puro

Vector name	pLenti-CMV-MCS-EF1a-mCherry-Puro
Cat No	VGB-V022-10C1
Size	8147 bp
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	mCherry

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
Note	For research use only, not for use in diagnostic procedure.

