


pLenti-mCherry-Puro+GPC3(h) plasmid

Cat No.: VGP-2719-33E2

Gene/insert name	Homo sapiens glypican 3
Official Symbol	GPC3
Alt name	DGSX;GTR2-2;MXR7;OCI-5;SDYS;SGB;SGBS;SGBS1
Insert size	1743 bp
Gene ID	2719
Accessions	NM_001164617.2
Tag	W/O
Backbone size	8138 bp
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	mCherry

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

