


pLenti-GFP-Puro+shRNA·NUP85(h)五靶点质粒套餐

Cat No.: VGP-79902-SHK2B

Gene/insert name	Homo sapiens nucleoporin 85
Official Symbol	NUP85
Alt name	FROUNT; NPHS17; Nup75
Gene ID	79902
Tag	W/O
Cloning site 5'	Age I
Cloning site 3'	EcoR I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	GFP
Promoter	U6
Loop	CTCGAG

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

