


pLenti-Puro+myc·STUB1(h) plasmid

Cat No.: VGP-10273-34J9

Gene/insert name	Homo sapiens STIP1 homology and U-box containing protein 1
Official Symbol	STUB1
Alt name	CHIP; HSPABP2; NY-CO-7; SCA48; SCAR16; SDCCAG7; UBOX1
Insert size	909 bp
Gene ID	10273
Accessions	NM_005861.4
Tag	myc
Backbone size	7409 bp
Cloning site 5'	Xho I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	W/O

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

