


pLenti-SV40T-T2A-Puro plasmid

Cat No.: VGP-B103-50F4

Gene/insert name	SV40
Official Symbol	SV40
Alt name	Large T antigen
Insert size	2127 bp
Gene ID	29031019
Accessions	YP_003708382.1
Tag	W/O
Backbone size	7095 bp
Cloning site 5'	EcoR I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	W/O

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

