


pLenti-Puro+Cd63(m)·NanoLuc plasmid

Cat No.: VGP-12512-35J8

Gene/insert name	Mus musculus CD63 antigen
Official Symbol	Cd63
Alt name	ME491; Tspan30
Insert size	714 bp
Gene ID	12512
Accessions	NM_001042580.2
Tag	NanoLuc
Backbone size	7361 bp
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	W/O

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

