


pLenti-Puro+HLA-B4601 plasmid

Cat No.: VGP-HLA-37A4

Gene/insert name	HLA-B4601
Official Symbol	HLA-B4601
Alt name	N/A
Insert size	1089 bp
Gene ID	W/O
Accessions	W/O
Tag	W/O
Backbone size	7508 bp
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	W/O

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

