


pLenti-Puro+GFP·APEX2·Flag plasmid

Cat No.: VGP-27301-33B4

Gene/insert name	APEX2
Official Symbol	APEX2
Alt name	APE2; APEXL2; XTH2; ZGRF2
Insert size	747 bp
Gene ID	27301
Accessions	NM_014481.4
Tag	Flag
Backbone size	8174 bp
Cloning site 5'	Xba I
Cloning site 3'	Xho I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	GFP

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

