


pLenti-mCherry-Puro+Myh9 plasmid

Cat No.: VGP-4627-33C9

Gene/insert name	Myh9
Official Symbol	MYH9
Alt name	BDPLT6; DFNA17; EPSTS; FTNS; MATINS; MHA
Insert size	2943 bp
Gene ID	4627
Accessions	NM_002473.6
Tag	W/O
Backbone size	8064 bp
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	mCherry

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

