


pLenti-eGFP-NeoR/KanR+RASEF(G892A,h) plasmid

Cat No.: VGP-158158-75A7

Gene/insert name	Homo sapiens RAS and EF-hand domain containing, mRNA
Official Symbol	RASEF(h)
Alt name	RAB45; TSG
Insert size	2223 bp
Gene ID	158158
Accessions	NM_152573.4
Tag	N/A
Backbone size	8986 bp
Cloning site 5'	EcoR I
Cloning site 3'	Not I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	NeoR/KanR
Fluorescent protein	eGFP

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

