


pLenti-Puro+shRNA-SRGN(h) 五靶点干扰质粒套餐

Cat No.: VGP-5552-SHK2A

Target gene name	Homo sapiens serglycin (SRGN), transcript variant 1, mRNA
Alt name	PPG; PRG; PRG1
Gene ID	5552
Plasmid size	7168 bp
Cloning site 5'	Age I
Cloning site 3'	EcoR I
Promoter	U6
Loop	CTCGAG
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	W/O

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.



产品信息

产品编号	产品名称	规格
VGW-5552-111J7	111J7_pLenti-Puro+shRNA·SRGN(h)-1	1E8 TU,1 mL
VGW-5552-111J8	111J8_pLenti-Puro+shRNA·SRGN(h)-2	1E8 TU,1 mL
VGW-5552-111J9	111J9_pLenti-Puro+shRNA·SRGN(h)-3	1E8 TU,1 mL
VGW-5552-112A1	112A1_pLenti-Puro+shRNA·SRGN(h)-4	1E8 TU,1 mL
VGW-5552-112A2	112A2_pLenti-Puro+shRNA·SRGN(h)-5	1E8 TU,1 mL
VGW-V107-23C9	23C9_pLenti-Puro+shRNA·NC	1E8 TU,1 mL

