


pLenti -Puro+ACVRL1(h) plasmid

Cat No.: VGP-0094-32D7

Gene/insert name	Homo sapiens activin A receptor like type 1
Official Symbol	ACVRL1
Alt name	ACVRLK1; ALK-1; ALK1; HHT; HHT2; ORW2; SKR3; TSR-I
Insert size	1512 bp
Gene ID	94
Accessions	NM_000020.3
Tag	N/A
Backbone size	7350 bp
Cloning site 5'	Xha I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	W/O

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

