


pLenti -RFP-Puro+NAT2(h) • His plasmid

Cat No.: VGP-10-001P1

Gene/insert name	Homo sapiens N-acetyltransferase 2, mRNA
Official Symbol	NAT2(h)
Alt name	AAC2; NAT-2; PNAT
Insert size	870 bp
Gene ID	10
Accessions	NM_000015.3
Tag	His
Backbone size	8115bp
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	RFP

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

