


pLenti - Puro+ MLPH (Sus scrofa) plasmid

Cat No.: VGP-100049680-51H9

Gene/insert name	Sus scrofa melanophilin (MLPH), mRNA
Official Symbol	MLPH (Sus scrofa)
Alt name	MLPH
Insert size	1734 bp
Gene ID	100049680
Accessions	NM_001098594.1
Tag	N/A
Backbone size	7362 bp
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	W/O

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

