


pLenti -Puro+ACTR3(h,K42A) plasmid

Cat No.: VGP-10096-001P2

Gene/insert name	Homo sapiens actin related protein 3, transcript variant 1, mRNA
Official Symbol	ACTR3(h)
Alt name	ARP3
Insert size	1257 bp
Gene ID	10096
Accessions	NM_005721.5
Tag	N/A
Backbone size	7362 bp
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	W/O

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

