


pLvX-Puro+TRAIP(h) • myc plasmid

Cat No.: VGP-10293-36C3

Gene/insert name	Homo sapiens TRAF interacting protein, mRNA
Official Symbol	TRAIP
Alt name	RNF206; SCKL9; TRIP
Insert size	1407 bp
Gene ID	10293
Accessions	NM_005879.2
Tag	myc
Backbone size	8099 bp
Cloning site 5'	Xho I
Cloning site 3'	Xba I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	W/O

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

