


pLenti-eGFP-NeoR+ Progerin(h) plasmid

Cat No.: VGP-0000-75B1

Gene/insert name	Homo sapiens progerin mRNA, complete cds
Official Symbol	Progerin
Alt name	N/A
Insert size	1845 bp
Gene ID	0000
Accessions	AY357727.1
Tag	N/A
Backbone size	9001 bp
Cloning site 5'	Xho I
Cloning site 3'	Not I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	eGFP

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

