


pLenti-GFP-Puro+Tigit(m) • Flag plasmid

Cat No.: VGP-100043314-001P1

Gene/insert name	Mus musculus T cell immunoreceptor with Ig and ITIM domains
Official Symbol	Tigit
Alt name	Vstm3
Insert size	726 bp
Gene ID	100043314
Accessions	NM_001146325.1
Tag	Flag
Backbone size	8239 bp
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	GFP

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

