


pLenti-GFP+ Isg15(m) • 3xFlag plasmid

Cat No.: VGP-100038882-001P1

Gene/insert name	Mus musculus ISG15 ubiquitin-like modifier (Isg15), mRNA
Official Symbol	Isg15
Alt name	100038882; G1p2; IGI15; IP17; Irfp; UCRP
Insert size	486 bp
Gene ID	100038882
Accessions	NM_015783.3
Tag	3xFlag
Backbone size	7590 bp
Cloning site 5'	Xba I
Cloning site 3'	Not I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	GFP

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

