


pLenti-GFP-Puro+LINC00894(h) plasmid

Cat No.: VGP-100272228-001P1

Gene/insert name	Homo sapiens EOLA2 divergent transcript
Official Symbol	LINC00894(h)
Alt name	EOLA2-DT
Insert size	3182 bp
Gene ID	100272228
Accessions	NR_027456.1
Tag	N/A
Backbone size	8158 bp
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	GFP

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

