


pLenti-GFP-Puro+ POU5F1(h) plasmid

Cat No.: VGP-100302349-001P1

Gene/insert name	Ovis aries POU class 5 homeobox 1, transcript variant X1, mRNA
Official Symbol	POU5F1(h)
Alt name	N/A
Insert size	1083 bp
Gene ID	100302349
Accessions	XM_004018968.5
Tag	N/A
Backbone size	8177 bp
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	GFP

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

