


pLenti-GFP-Puro+ MKLN1-AS(h) plasmid

Cat No.: VGP-100506881-001P1

Gene/insert name	Homo sapiens MKLN1 antisense RNA, long non-coding RNA
Official Symbol	MKLN1-AS(h)
Alt name	MKLN1-AS2
Insert size	2046 bp
Gene ID	100506881
Accessions	NR_125364.1
Tag	N/A
Backbone size	8177 bp
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	GFP

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

