


pLenti-ffluc2-Puro+IL10(m) • 3xFlag plasmid

Cat No.: VGP-16153-43F9

Gene/insert name	Mus musculus interleukin 10 (Il10), mRNA
Official Symbol	IL10
Alt name	CSIF; If2a; Il-10
Insert size	537 bp
Gene ID	16153
Accessions	NM_010548.2
Tag	3xFlag
Backbone size	7987 bp
Cloning site 5'	EcoR I
Cloning site 3'	Not I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	W/O

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

