


pLenti-GFP-Puro+IL10(m) • His8 plasmid

Cat No.: VGP-16153-58D6

Gene/insert name	Mus musculus interleukin 10 (Il10), mRNA
Official Symbol	IL10
Alt name	CSIF; If2a; Il-10
Insert size	537 bp
Gene ID	16153
Accessions	NM_010548.2
Tag	His8
Backbone size	7987 bp
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	GFP

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

