


pLenti-BSD+shRNA·Chd1l(m)五靶点干扰质粒套餐

Cat No.: VGP-68058-SHK2D

Target gene name	Mus musculus chromodomain helicase DNA binding protein 1-like
Alt name	4432404A22Rik; Alc1; Snf2p
Gene ID	68058
Plasmid size	6991 bp
Cloning site 5'	Age I
Cloning site 3'	EcoR I
Promoter	U6
Loop	CTCGAG
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Blasticidin
Fluorescent protein	W/O

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

