


pLenti-Puro+myc • mTrkA CAR PART A+B plasmid

Cat No.: VGP-0000-12B9

Gene/insert name	mTrkA CAR PART A+B
Official Symbol	N/A
Alt name	N/A
Insert size	N/A
Gene ID	N/A
Accessions	N/A
Tag	myc
Backbone size	N/A
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	eGFP+mCherry

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

