


pLenti-GFP-Puro+shRNA·COL1A2(h) 五靶点干扰质粒套餐

Cat No.: VGP-1278-SHK2B

Target gene name	Homo sapiens collagen type I alpha 2 chain, mRNA
Alt name	EDSARTH2; EDSCV; OI4
Gene ID	1278
Plasmid size	7970 bp
Cloning site 5'	Age I
Cloning site 3'	EcoR I
Promoter	U6
Loop	CTCGAG
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	GFP

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

