


## pLenti-KanR+shRNA·Med19(m) 三靶点干扰质粒套餐

**Cat No.:** VGP-381379-Med19

<b>Target gene name</b>	Mus musculus mediator complex subunit 19 (Med19)
<b>Alt name</b>	2410018M14Rik; 3110040A13Rik; LCMR1
<b>Gene ID</b>	381379
<b>Plasmid size</b>	7632bp
<b>Cloning site 5'</b>	Age I
<b>Cloning site 3'</b>	EcoR I
<b>Promoter</b>	U6
<b>Loop</b>	CTCGAG
<b>Growth strain(s):</b>	Stbl 3
<b>High or low copy</b>	High
<b>Bacterial resistance(s)</b>	Ampicillin
<b>Growth temperature</b>	37 °C
<b>Vector type</b>	Lentiviral, Mammalian Expression
<b>Selectable markers</b>	Puromycin
<b>Fluorescent protein</b>	W/O

### Map



<b>Buffer</b>	TE buffer
<b>Storage</b>	Store at -20°C up to 24 months
<b>The plasmid is ready for</b>	Transfection
<b>Note</b>	For research use only, not for use in diagnostic procedure.

产品信息



产品编号	产品名称	规格
VGX-381379-113B7	pLenti-Puro+shRNA • Med19(m)-1	1E8 TU,1 mL
VGX-381379-113B8	pLenti-Puro+shRNA • Med19(m)-2	1E8 TU,1 mL
VGX-381379-113B9	pLenti-Puro+shRNA • Med19(m)-3	1E8 TU,1 mL
VGX-381379-113C1	pLenti-Puro+shRNA • Med19(m)-4	1E8 TU,1 mL
VGX-381379-113C2	pLenti-Puro+shRNA • Med19(m)-5	1E8 TU,1 mL
VGP-V000-48C3	pLenti-Puro+shRNA • MCS(U6)载体	1E8 TU,1 mL

